

# Correlations between Cefiderocol Broth Microdilution MICs and Disk Diffusion Inhibitory Zone Diameters among Target Gram-Negative Organisms

Masakatsu Tsuji <sup>1)</sup>, Meredith Hackel <sup>2)</sup>, Naoki Koira <sup>1)</sup>, Yoshinori Yamano <sup>1)</sup>, Roger Echols <sup>3)</sup> and Daniel F Sahn <sup>2)</sup>  
1) Shionogi & Co., Ltd., Osaka JP, 2) IHMA, Inc. Schaumburg, IL, US, 3) ID3C, LLC, Easton CT, US

Contact information:  
Masakatsu Tsuji  
Address: 3-1-1, Futaba-cho, Toyonaka, Osaka, 561-0825, JAPAN  
TEL: +81-6-6331-8272 FAX: +81-6-6331-8612  
E-mail address: masakatsu.tsuji@shionogi.co.jp

## ABSTRACT

**Background:** Cefiderocol (S-649266) is a novel siderophore cephalosporin with potent activity against a wide variety of Gram-negative bacteria, including carbapenem-resistant strains. In order to evaluate the utility of disk diffusion for in vitro susceptibility testing of cefiderocol, this study was conducted to determine the correlation between broth microdilution MICs and disk diffusion inhibitory zone diameters for this drug.

**Materials and Methods:** In all, 1319 Gram-negative clinical isolates were included in this study (807 Enterobacteriaceae, 173 *Pseudomonas aeruginosa*, 254 *Acinetobacter baumannii*, 3 *Burkholderia cepacia*, and 82 *Stenotrophomonas maltophilia*). The study isolates represented a global collection obtained during 2014 and 2015. All isolate identifications were confirmed and susceptibility testing was performed central-ly (IHMA, Inc., USA). Susceptibility testings were conducted according to the Clinical and Laboratory Standard Institute guidance. For broth microdilution, iron-depleted CAMHB (ID-CAMHB) was used. Disk diffusion tests using the disk containing 30 µg cefiderocol were performed on unsupplemented standard Mueller-Hinton agar (MHA). In addition, the potential effect of the addition of iron into MHA for disk diffusion testing was evaluated with two QC strains (*E. coli* ATCC 25922 and *P. aeruginosa* ATCC 27853).

**Results:** For the two QC strains (*E. coli* ATCC 25922 and *P. aeruginosa* ATCC 27853) the addition of iron at 1 mg/L to MHA did not notably impact the zone diameters obtained. Because the iron content in MHA is approximately 0.3 mg/L, it was determined that disk diffusion studies could be conducted using standard MHA. Correlation between MIC values in ID-CAMHB and disk zones was relatively high for all test isolates. The R<sup>2</sup> values observed for Enterobacteriaceae, *A. baumannii*, *P. aeruginosa*, and *S. maltophilia* were 0.70, 0.81, 0.52, and 0.65, respectively. The relatively low R<sup>2</sup> values observed for *P. aeruginosa* and *S. maltophil-ia* may not reflect true correlation due to the lack of resistant strains.

**Conclusions:** The disk diffusion studies using a 30 µg of cefiderocol disk and standard MHA showed that the zone diameters obtained in these studies have a good correlation with the MIC determined using ID-CAMHB. Disk diffusion test was validated as an alternative method for testing cefiderocol MIC for Gram-negative bacteria.

## INTRODUCTION

Cefiderocol is a novel parenteral catechol-substituted siderophore cephalosporin that is active against carbapenem-resistant Gram-negative bacteria (1, 2, 3). This study was performed to determine the effect of iron on the *in vitro* antibacterial activity and for correlation of cefiderocol MIC and disk diffusion method.

## MATERIALS AND METHODS

### Test organisms

In all 1319 Gram-negative bacteria isolates (807 Enterobacteriaceae, 254 *Acinetobacter baumannii*, 173 *Pseudomonas aeruginosa*, 82 *Stenotrophomonas maltophilia*, and 3 *Burkholderia cepacia*) were analyzed. These clinical isolates were obtained from IHMA, Inc. Two ATCC quality control (QC) organisms were tested: *Escherichia coli* ATCC2522 and *P. aeruginosa* ATCC27853.

### Compound

Cefiderocol and cefepime were used. Cefiderocol disks were manufactured by MAST (Merseyside, UK) and Bio-Rad (Marnes-la-Coquette,FR).

### MIC determination

The appropriate medium for the disk diffusion study of cefiderocol has been approved by Clinical and Laboratories Standard Institute (CLSI) guidelines to use standard MHA. A standardized inoculum of bacteria is swabbed onto the surface of a Mueller-Hinton agar plate. Filter paper disks with cefiderocol are placed on the agar. After overnight incubation in a 35-37°C ambient incubator for 18-20 hours, the diameter of the zone of inhibition is measured around each disk. CLSI approved cefiderocol disk QC range was used as follows: *E. coli* ATCC2522 (25-31 mm) and *P. aeruginosa* ATCC27853 (22-31 mm). To check the correlation of both MIC determination methods, broth microdilution testing was done by using iron-depleted cation-adjusted Muller-Hinton Broth (ID-CAMHB) according to CLSI (4, 5).

## RESULTS

Table 1. List of tested isolates

| Species                             | N           | Percent of total (%) |
|-------------------------------------|-------------|----------------------|
| <b>Enterobacteriaceae</b>           | <b>807</b>  | <b>61.2</b>          |
| <i>Citrobacter freundii</i>         | 26          | 2.0                  |
| <i>Citrobacter koseri</i>           | 14          | 1.1                  |
| <i>Enterobacter aerogenes</i>       | 44          | 3.3                  |
| <i>Enterobacter asburiae</i>        | 17          | 1.3                  |
| <i>Enterobacter cloacae</i>         | 149         | 11.3                 |
| <i>Escherichia coli</i>             | 167         | 12.7                 |
| <i>Klebsiella oxytoca</i>           | 45          | 3.4                  |
| <i>Klebsiella pneumoniae</i>        | 250         | 19.0                 |
| <i>Proteus mirabilis</i>            | 2           | 0.2                  |
| <i>Serratia marcescens</i>          | 93          | 7.1                  |
| <b>Non-Enterobacteriaceae</b>       | <b>512</b>  | <b>38.8</b>          |
| <i>Acinetobacter baumannii</i>      | 254         | 19.3                 |
| <i>Burkholderia cepacia</i>         | 3           | 0.2                  |
| <i>Pseudomonas aeruginosa</i>       | 173         | 13.1                 |
| <i>Stenotrophomonas maltophilia</i> | 82          | 6.2                  |
| <b>Total</b>                        | <b>1319</b> | <b>100</b>           |


Table 2. Effect of iron on cefiderocol disk zone diameter against Gram-negative bacteria

| Test Medium                        | Diameters (mm)            |          |                                 |          |
|------------------------------------|---------------------------|----------|---------------------------------|----------|
|                                    | <i>E. coli</i> ATCC 25922 |          | <i>P. aeruginosa</i> ATCC 27853 |          |
|                                    | Cefiderocol               | Cefepime | Cefiderocol                     | Cefepime |
| MHA Only                           | 30 – 31                   | 37       | 29 – 30                         | 30 – 31  |
| MHA + FeCl <sub>3</sub> (10 mg/L)  | 27 – 29                   | 37       | 18 – 19                         | 31       |
| MHA + FeCl <sub>3</sub> (3 mg/L)   | 28 – 29                   | 37       | 21 – 22                         | 30 – 31  |
| MHA + FeCl <sub>3</sub> (1 mg/L)   | 30                        | 37 – 38  | 26 – 28                         | 30 –31   |
| MHA + FeCl <sub>3</sub> (0.3 mg/L) | 30                        | 37       | 28 – 30                         | 31       |

Cefiderocol disk: 30 µg/disk (Bio-rad). CFPM disk 30 µg/disk (Becton, Dickinson and Company).  
Medium: Difco MHA, BBL MHA (Becton, Dickinson and Company)

Little effect of iron concentration on disk diameter of cefiderocol was observed even when 10 mg/L of iron for *E. coli* ATCC 25922 and 1 mg/L of iron for *P. aeruginosa* ATCC 27853 were added to the MHA.

Figure 1. Cefiderocol MIC values compared to cefiderocol disk zone of inhibition (30 µg/disk) against 1,319 Gram-negative bacteria



There is good correlation between disk diffusion and BMD method for cefiderocol MIC when multiple bacterial species (a total of 1,319 strains including 807 Enterobacteriaceae, 254 *A. baumannii*, 173 *P. aeruginosa*, 82 *S. maltophilia*, 3 *B. cepacia*) were used.

Figure 2. Cefiderocol MIC values compared to cefiderocol disk zone of inhibition (30 µg/disk) against 807 Enterobacteriaceae




Figure 3. Cefiderocol MIC values compared to cefiderocol disk zone of inhibition (30 µg/disk) against 173 *P. aeruginosa*




Figure 4. Cefiderocol MIC values compared to cefiderocol disk zone of inhibition (30 µg/disk) against 254 *A. baumannii*




Figure 5. Cefiderocol MIC values compared to cefiderocol disk zone of inhibition (30 µg/disk) against 82 *S. maltophilia*




Table 3. Correlation between cefiderocol MIC values in BMD method and disk zones diameter.

| Organism group                      | R2 values |
|-------------------------------------|-----------|
| Enterobacteriaceae                  | 0.70      |
| <i>Acinetobacter baumannii</i>      | 0.81      |
| <i>Pseudomonas aeruginosa</i>       | 0.52      |
| <i>Stenotrophomonas maltophilia</i> | 0.65      |

• Correlation between MIC values in ID-CAMHB and disk zones was high for Enterobacteriaceae and *A. baumannii*, with R<sup>2</sup> values of 0.70 and 0.81, respectively, indicating that a disk test is very feasible for use with this organism group.  
• R<sup>2</sup> values were reduced for *P. aeruginosa* (R<sup>2</sup>=0.52) and *S. maltophilia* (R<sup>2</sup>=0.65). These lower values may not reflect true correlation due to the lack of resistant strains.

## CONCLUSIONS

The disk diffusion studies using a 30 µg of cefiderocol disk and standard MHA showed that the zone diameters obtained in these studies have a good correlation with the MIC determined using ID-CAMHB.

## REFERENCES

- M A. Hackel, M. Tsuji, Y. Yamano, R. Echols, J A. Karlowsky, and D F Sahn. *In Vitro* Activity of the Siderophore Cephalosporin, Cefiderocol, against Carbapenem-Nonsusceptible and Multidrug-Resistant Isolates of Gram-Negative Bacilli Collected Worldwide in 2014 to 2016: AAC 2016 February 62(2): e01968-17.
- M A. Hackel, M. Tsuji, Y. Yamano, R. Echols, J A. Karlowsky, and D F Sahn. In Vitro Activity of the Siderophore Cephalosporin, Cefiderocol, against Carbapenem-Nonsusceptible and Multidrug-Resistant Isolates of Gram-Negative Bacilli Collected Worldwide in 2014 to 2016: AAC 2018 February 62(2): e01968-17.
- A. Ito, T. Sato, M. Ota et al. *In Vitro* Antibacterial Properties of Cefiderocol, a Novel Siderophore Cephalosporin, against Gram-Negative Bacteria AAC 2018 February 62(2): e014-17.
- Clinical Laboratory Standards Institute (CLSI), 2015. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standards -- Tenth Edition. CLSI document M07-A10 (ISBN 1-56238-988-2). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA.
- Clinical and Laboratory Standards Institute (CLSI), 2016. Performance Standards for Antimicrobial Susceptibility Testing; Twenty-Fifth Informational Supplement. CLSI Document M100-S26 (ISBN 1-56238-923-8). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA.